AQUOLAB[®]

AQUOLAB[®] BASIC

- IT Manuale di istruzioni
- EN User's manual
- ES Manual de instrucciones
- P Manual de instruções

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AQUOLAB® BASIC - DESCRIPTION OF COMPONENTS

- a Nozzles (4 per device included)
- b Nozzle container (box separate from the dental water jet)
- c ON/OFF button
- d Ozone activation and regulation button
- e Water pressure selection button
- f Push-button to turn the water flow ON/OFF
- g Device handle
- h Tank (300 ml)
- i Tank cover
- j Reset button


FIGURE 1

Aquolab Basic is a certified electro-medical device. Read the warnings and instructions for use carefully.

Any other use than the intended one (oral hygiene) constitutes an improper use. The manufacturer cannot be held liable for any damage caused by improper, incorrect and/or unreasonable use or if the device is connected to electrical systems which do not comply with safety regulations in force.

Aquolab is a device for the whole family, but children should only use it under adult supervision. It should not be used by children under 6 years of age.

The Aquolab dental water jet mixes water and ozone, directly produced inside the device in controlled quantities well below the European standard safety levels, in order to provide you with better dental hygiene and oral care. Thanks to its antibacterial, antimicrobial, immunomodulating and anti-inflammatory properties, ozone is able to rapidly decompose bacteria and viruses, preventing inflammatory processes and promoting healing. Aquolab's jet of water and ozone eliminates bacteria and viruses in just a few seconds' time, gently, painlessly and without leaving any residues, because the ozone transforms into oxygen after a few seconds. Aquolab works effectively using regular potable water. It is recommended that you use Aquolab after normal tooth brushing. For proper oral hygiene, we recommend using it once a day, preferably in the evening.

ACTIVATION

You must read the instruction manual and safety rules carefully before use.

- Place Aquolab on a horizontal surface.
- Connect the power supply charger to an easily accessible outlet.
- Connect to an electrical system constructed in accordance with the safety standards.
- Now your Aquolab Basic is ready.

ACCESSORIES

As standard equipment you will receive 4 nozzles: 3 with a 0.6 mm diameter (coloured nozzles) and 1 with a 0.8 mm diameter (white nozzle). The package you will receive also contains a blister package with 3 vials of de-scaling solution for ordinary maintenance (see the section on Maintenance).

Nozzles for periodic replacement, with 0.6 mm diameter and 0.8 mm diameter, as well as the de-scaling solution, are available on the website www.aquolab.com.

In the package you will also find the handle, already correctly inserted into the device. The handle is detachable. If it is damaged or requires replacement, you can replace it without having to replace the entire device.

HOW TO USE THE AQUOLAB® DENTAL WATER JET

1. When using for the first time or after a long period of inactivity, perform a vacuum cycle before using your dental water jet
2. The device switches off automatically after 1 min/1.5 min to prevent water-free operation
3. Fill the tank with water, about 300 ml.
4. Insert your personal nozzle into the handle.
5. When you push the ON button (c), the LED will flash to remind you to make sure that there is water in the tank (h). If there is water in the tank, proceed to step 6.
6. To turn the water flow ON/OFF, press the dedicated button. (e)
7. To turn the ozone flow ON/OFF, press the dedicated button. (d)
8. Press the button activating the water flow. (f)
9. For children under 10 years of age the 0.8 mm diameter nozzle is strongly recommended.
10. If you have high dentine sensitivity, use the 0.8 mm diameter nozzle.

HOW TO REGULATE WATER AND OZONE LEVELS

Once you have switched on your device, you can regulate/modify the levels of water and ozone by following these steps:

- Press and hold down the C and D buttons SIMULTANEOUSLY for 3 seconds.
- The orange light (J) will come on.
- To set the water level, press button C. The LED will flash once or twice (level 1 or level 2). To change the level, press button C for a second time.
- To set the ozone level, press button D. The LED will flash once or twice (level 1 or level 2). To change the level, press button D for a second time.
- Once you have set the desired levels, unplug the device, wait for a few seconds, then plug it back in. Now your Aquolab will operate and it will save the levels you have selected.

These are just general recommendations. For specific needs, contact our Customer Service by calling number: +39 02 366 86 416 or writing an email to: assistenza@aquolab.com

ALARM NOTIFICATION AND RESTORATION

If a fault occurs during operation, the device stops the flow and signals the fault condition with a flashing LED alarm (j).

To restore normal operation, press the ON/OFF button (c). If the fault persists, turn the device off and then restart it. If the fault signal still persists, then the device requires specialist assistance (contact Customer Service).

IMPORTANT

Periodically check the condition of the external and visible parts of Aquolab. If you notice that the housing, the charger or the cable is damaged, do not use the device. Contact customer service (+39 02 366 86 416 – assistenza@aquolab.com).

Do not use the charger or cables if they are damaged. Replacements are available for all of the external parts. Aquolab will automatically warn you of any malfunctions that occur in its internal components. If the solution to the problem cannot be found in this manual, contact customer service.

IMPORTANT WARNINGS

- Use only with the provided EM1024RE-24 power supply charger.
- Fill the tank with water. DO NOT use the device without water or if it fails to dispense water.
- Use fresh water (not stagnant) every time.
- Use water at room temperature only (avoid warm or cooled water).
- Use a different nozzle for each user.
- Its use is intended for the oral cavity only.
- DO NOT use the device if its external parts or the charger are damaged.
- Use only for the intended use.
- Do not modify the device in any way.
- Pay particular attention to use by children to prevent possible strangulation with the power cord or handle tube
- De-scale the device regularly, using only the specific product indicated.

DIRECTIONS FOR ELECTROMAGNETIC COMPATIBILITY

Since Aquolab is an electro-medical product, you need pay attention to the following tables, in order to let the device operate safely in accordance with the information about electro-medical machines.

It is important to remember that portable and mobile radio communication can affect the device's operation. It is recommended that any electro-medical device – therefore, including Aquolab – not be used near or stacked with other equipment.

TABLE 1 - GUIDELINES AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS

RF emissions CISPR 11	Group 1	The device "AQUOLAB BASIC" uses RF energy only for its own internal operation. Therefore, its RF emissions are very low and are not likely to cause any interference to nearby electronic equipment.
RF emissions CISPR 11	Class B	
Emission of harmonics IEC 61000-3-2	Class A	
Emission of voltage fluctuations and flicker IEC 61000-3-3	Compliance	

TABLE 2 - GUIDELINES AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC INTERFERENCE RESISTANCE (IMMUNITY)

Electrostatic Discharge (ESD) IEC 61000-4-2	± 6kV contact ± 8 kV air	± 6kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Fast transient/electrical bursts IEC 61000-4-4	± 2kV for power supply lines	± 2kV for power supply lines	The mains power supply quality should be that of a typical commercial or hospital environment.
Surges IEC 61000-4-5	± 1 kV differential mode	± 1 kV differential mode	The mains power supply quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations of the power supply lines IEC 61000-4-11	< 5% U_i (> 95% dip in U_i) for ½ period 40% U_i (60% dip in U_i) for 5 periods 70% U_i (30% dip in U_i) for 25 periods < 5% U_i (> 95% dip in U_i) for 5s	< 5% U_i (> 95% dip in U_i) for ½ period 40% U_i (60% dip in U_i) for 5 periods 70% U_i (30% dip in U_i) for 25 periods < 5% U_i (> 95% dip in U_i) for 5s	The mains power supply quality should be that of a typical commercial or hospital environment. If the user of the device "AQUOLAB BASIC" requires continuous operation even in the case of power supply interruptions, it is recommended to connect the device "AQUOLAB BASIC" to an uninterruptible power supply or a battery.
Power frequency (50 Hz/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	-

TABLE 4 - RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE
RF TELECOMMUNICATIONS DEVICES AND THE "AQUOLAB BASIC"

Maximum output power rating of transmitter	Separation distance according to frequency of transmitter m		
	From 150 kHz to 80 MHz $d = 1.2 P$	From 80 MHz to 800 MHz $d = 1.2 P$	From 800 MHz to 2.5 GHz $d = 2.3 P$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

TABLE 3 - GUIDELINES AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC INTERFERENCE RESISTANCE
(08-04-00172)

Interference Resistance (Immunity) Test	IEC 60601 EN 60601-1-2	Compliance Level	Electromagnetic Environment – Guidelines
Conducted RF EN 61000-4-6	3 V eff From 150 kHz to 80 MHz	3 V eff From 150 kHz to 80 MHz	Portable and mobile RF communications devices should be used no closer to "AQUOLAB BASIC" or any of its components (including cables) than the recommended separation distance determined using the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2\sqrt{P}$ 150 kHz to 80 MHz Recommended separation distance: $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz Recommended separation distance: $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the information provided by the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, must be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:
Radiated RF EN 61000-4-3	3 V/m From 80 MHz to 2.5 GHz	3 V/m From 80 MHz to 2.5 GHz	

AQUOLAB® BASIC MAINTENANCE

Disconnect the device from the supply mains before any maintenance activity

ORDINARY CLEANING

Dry the water tank (h) with a cloth after every use to prevent limestone scale build-up.

ORDINARY REPLACEMENTS

It is recommended that the nozzles be replaced every 6-8 months.

ORDINARY MAINTENANCE – ANTI-SCALING TREATMENT

Aquolab Basic alerts you to the need for ordinary maintenance through the simultaneous flashing of the three LEDs. Remember that, at least once a month or when Aquolab gives you this signal, you need to run a de-scaling cycle with the specific de-scaling solution supplied in vials (Aquolab De-scaling solution) which you can purchase directly from the www.aquolab.com website. The de-scaling cycle allows you to remove scale build-up resulting from calcium and magnesium salts in water. This will maintain continuing device efficiency. To avoid damaging the device, do not use products other than the Aquolab de-scaling kit.

Follow these steps to correctly use the Aquolab De-scaling solution:

- Shake the vial before using it.
- Pour the contents into the dental water jet tank and fill it with water.
- Switch the device on and let it dispense water for approximately 15 sec.
- Wait 5 minutes and repeat this operation for a further 15 sec.
- Repeat the operation until the tank is completely empty.
- Rinse the tank by filling it with water only and letting it flow out at regular intervals until the tank is completely empty.

To reset the "Ordinary Maintenance" alarm and make Aquolab Basic ready for use, you need to press the reset button again (see Figure 1).

HOW TO REPLACE A DAMAGED HANDLE

Carefully detach the cable on the side of the dental water jet so as not to damage the outlet fitting.

Replace the handle with a new one, carefully inserting the cable into the outlet fitting.

EXTRAORDINARY MAINTENANCE – REPLACEMENT OF THE OZONE PRODUCER

The production of ozone lasts approximately 3 years (it may vary depending on number of use cycles). In this case, the device will alert you to the need for “Extraordinary Maintenance” through the sequential flashing of the LEDs. The device must be sent back to Customer Service so that the Ozone Producer can be replaced. It is not possible to independently reset the “Extraordinary Maintenance” alarm and the device will be able to work again only after the Ozone Producer has been replaced.

IMPORTANT

- Always remember to empty the hydraulic circuit before any transport or long-term storage of the device.
- If the device is going to be returned for extraordinary maintenance, remove the interchangeable nozzles before returning.
- When requesting a replacement handle, be very careful not to damage the outlet fitting when removing the old handle.
- **These are fundamental requirements. Non-compliant returns may be refused.**

DISPOSAL

Pay particular attention to the regulations on waste sorting and disposal currently in force. The water jet device must not be thrown away as household waste, but has to be handed over at authorized municipal landfills.

The external and replaceable parts of the water jet, such as nozzles and handle, can be recycled normally.

TECHNICAL SPECIFICATIONS

Manufacturer	EB2C srl - Via Savona, 94 - 20144 Milan - Italy
Model	Aquolab Basic
Voltage supply	100-240 V 50/60 Hz 36 W
Class	II
Applied part	Handle (type B)
Weight	0.99 Kg
Dimensions	190 mm x 140 mm x h150 mm
Usage temperature	+5°C ÷ +40°C
Intermittent operation	3 min ON / 1 min OFF
Transport temperature	-25°C ÷ +50°C

EXPLANATION OF LABELLING SYMBOLS

	Type B applied part
	It is obligatory to read this user's manual
	CE marking
	Selective waste sorting is obligatory
	Do not reuse (disposable device – just for some handles)
IP21	Protected against access with a finger to dangerous voltage parts Protected against the entry of vertical water drops

WARRANTY CONDITIONS

- The device is covered by warranty, valid 2 years from the date of purchase, for any faults in workmanship or materials used.
- Under this warranty, Aquolab will replace and/or repair originally defective components free of charge.
- The warranty does not cover the accessories supplied with the device or parts subjected to normal wear.
- The device must be repaired at authorized technical service centres only.
- The device must be sent to the authorized service centre within 8 days after the fault is noticed.
- The user is responsible for shipping costs.
- The user shall bear the expense for any repairs not covered by the warranty conditions.
- The warranty will be rendered void if the device is misused or abused, if the required regular maintenance (anti-scaling treatment) has not been performed, or if the damage cannot be attributed to the manufacturer.
- The warranty does not imply any compensation for any kind of damages to people or property, either direct or indirect, that occurs when the device is not operating correctly.
- The warranty is valid as from the date of purchase, certified by the cash register receipt or the invoice.
- The lack of the properly filled in warranty coupon accompanied by the above receipt/invoice voids the warranty.

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This medical device complies with the provisions of directive 93/42/EEC and its subsequent amendments.

Italian Ministry of Health Registration No. 1573613

Aquolab® is a registered trademark
No. 1231275 – EU 13245535

Industrial Patent No. 0001420117

Duplication is prohibited. Partial reproduction, in any form, is also prohibited unless previously authorized in writing by the Manufacturer.

Technical specifications, description, standard accessories supplied and product appearance may change without prior notice.

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